

**easywrap
CUSTOM**

| | | |
|---------------|-------------------|--------------------|
| Patient Name | Order No. | Repeat Garment No. |
| Date Measured | Clinic / Hospital | |
| Measured by | Telephone | E-mail |

1 QUANTITY

| FOOT | Right | Left |
|-----------------|-------|------|
| Light EL-FC | | |
| Strong ES-FC | | |



| LEG | Right | Left |
|-----------------|-------|------|
| Light EL-LC | | |
| Strong ES-LC | | |

| KNEE ONLY | Right | Left |
|-----------------|-------|------|
| Light EL-KC | | |
| Strong ES-KC | | |

| THIGH ONLY | Right | Left |
|-----------------|-------|------|
| Light EL-TC | | |
| Strong ES-TC | | |

MEASUREMENTS

| | THIGH | KNEE | LEG | FOOT |
|------------------------|---|------|-----|------|
| g | Just below the groin | | | |
| f | Mid-thigh | | | |
| e¹ | Bottom of the thigh | | | |
| e¹-g | Follow the contours of the limb | | | |
| e¹ | Bottom of the thigh | | | |
| e | Over knee | | | |
| d | 2 finger width below knee crease | | | |
| d-e¹ | Follow the contour of the limb | | | |
| c | Widest point of calf | | | |
| b² | Mid-point between b ¹ and c | | | |
| b¹ | Point at which the achilles tendon changes into calf muscle | | | |
| b | Ankle | | | |
| b-d | Follow the contour of the limb | | | |
| h | Heel | | | |
| a¹ | Mid foot | | | |
| a | Base of Toes | | | |
| A-b | Straight distance | | | |
| a-heel | Straight distance | | | |



4 ADDITIONAL LINERS

| | | STANDARD c circumference < 60cm | EXTRA WIDE c circumference > 60cm |
|--|------------|------------------------------------|--------------------------------------|
| PAIR OF A-d LINERS INCLUDED WHEN YOU ORDER FOOT & LEG | A-d | E-L-AD | E-L-ADW |
| PAIR OF A-g LINERS INCLUDED WHEN YOU ORDER FOOT, LEG & THIGH | A-g | E-L-AG | E-L-AGW |

5 COMMENTS

CAUTION: Assess shape of limb and ensure sufficient padding is used to protect vulnerable areas. easywrap Strong fabric should not be used on ankles under 18cm. For ankle measurements under 18cm apply enough padding to bring the ankle measurement to above 18cm. For ankles between 18cm and 22cm easywrap Strong should only be applied after specialist referral and under supervision. For ankles under 18cm easywrap Light may be used after specialist referral and under supervision.