

HOW TO USE THIS FORM

- Attach to a prescription for 'HADDENHAM CUSTOM GARMENTS'
- Sections ①, ②, ③ & ⑥ are mandatory. Please fill in as appropriate.
- Sections ④, ⑤ and ① only require filling out if you require any of these options.
- Where Haddenham Custom Made garments are prescribed including specifically prescribed additional options, the dispenser must ensure that full details including order codes are clearly endorsed on the Drug Tariff. The additional options for the base garments do not attract any additional prescription charges.
- One prescription charge is payable where one or more of the same base garment is supplied. More than one prescription charge is payable if different base garments (Including different compression of the same garment type) are ordered.

NOTE: Where Haddenham Custom Made garments are prescribed including specifically prescribed additional options, the dispenser must ensure that full details including order codes are clearly endorsed on the Drug Tariff. The additional options for the base garments do not attract any additional prescription charges.

One prescription charge is payable where one or more of the same base garment is supplied. More than one prescription charge is payable if different base garments (Including different compression of the same garment type) are ordered (see part XVI sub-paragraph 12.13 of drug tariff).

RTFX a softer fabric for the upper limb

GOI DPIINKT a firmer fabric for the upper limb



Hadden Class

Class Class

NGERS with FINA

nham Custom PERTEX	CODE	Haddenham Custom PERTI
1 - 18-21 mmHg	PX-61	Class 1 - 18-21 mmHg
2 - 23-32 mmHg	PX-62	Class 2 - 23-32 mmHg
3 - 34-46 mmHg	PX-63	Class 3 - 34-46 mmHg

Class 2 - 23-32 mmHg	GP-6
Class 3 - 34-46 mmHg	GP-6



addenham Custom PERTEX	CODE	Haddenham Custom PER
lass 1 - 18-21 mmHg	PX-64	Class 1 - 18-21 mmHg
lass 2 - 23-32 mmHg	PX-65	Class 2 - 23-32 mmHg
lass 3 - 34-46 mmHg	PX-66	Class 3 - 34-46 mmHg

Haddenham Custom GOLDPUNK

GP-62	Class 2 - 23-32 mmHg	GP-65
GP-63	Class 3 - 34-46 mmHg	GP-66

Haddenham	Custom	GOLDPUNKT

Class 2 - 23-32 mmHg	GP-72
Class 3 - 34-46 mmHg	GP-73

PX-71

PX-73



Haddennam Custom PERTEX	CODE
Class 1 - 18-21 mmHg	PX-74
Class 2 - 23-32 mmHg	PX-75
Class 3 - 34-46 mmHg	PX-76

GLOVE with FINGE

Haddenham Custom PERTEX	CODE
Class 1 - 18-21 mmHg	PX-7
Class 2 - 23-32 mmHg	PX-7
Class 3 - 34-46 mmHg	PX-7

റ	PT	IN	NIC	
U	ГΙ	IU	IVO)

PERTEX and GOLDPUNKT	CODE
SILICONE GRIP TOP	CM-8
ZIP	CM-8
SHOULDER CAP	CM-8

Haddenham Custom GOLDPUNKT

Class 2 - 23-32 mmHg	GP-75
Class 3 - 34-46 mmHg	GP-76

Haddennam Custom GULDPUNK	ı
Class 2 - 23-32 mmHg	GP-78
Class 3 - 34-46 mmHg	GP-79

Responsibility for the selection and application of products rests with the professional staff treating the patient.

Haddenham Healthcare does not take any responsibility for contraindicated use of these garments or for inaccurate measurements.