

## HOW TO USE THIS FORM

- Attach to a prescription for 'HADDENHAM CUSTOM GARMENTS'
- Sections (1), (2), (3) & (4) are mandatory. Please fill in as appropriate. (N.B. If the garment is 'footless' you will not need to fill out section (3).)
- Sections (5), (6) and (7) only require filling out if you require any of these options.
- For invoicing purposes garments with measurements up to 'd' are BELOW KNEE. Garments with measurements above 'd' but without a waist attachments and/or panty measurements are THIGH HIGH. Garments with a 'waist attachment' are THIGH HIGH WITH WAIST ATTACHMENT. Garments with panty measurements ('t', 'j') are PANTYHOSE.

NOTE: Where Haddenham Custom Made garments are prescribed including specifically prescribed additional options, the dispenser must ensure that full details including order codes are clearly endorsed on the Drug Tariff. The additional options for the base garments do not attract any additional prescription charges.

One prescription charge is payable where one or more of the same base garment is supplied. More than one prescription charge is payable if different base garments (Including different compression of the same garment type) are ordered (see part XVI sub-paragraph 12.13 of drug tariff).

<b>PERTEX</b> a softer fabric for the lower limb <b>GOL</b>			PUNKT	a firmer fabric for the low	er limb
BELOW KNEE		THIGH HIGH			
Haddenham Custom PERTEX	CODE	Haddenham Custom PERTEX	CODE		
Class 1 - 18-21 mmHg	PX-21	Class 1 - 18-21 mmHg	PX-31		
Class 2 - 23-32 mmHg	PX-22	Class 2 - 23-32 mmHg	PX-32		
Class 3 - 34-46 mmHg	PX-23	Class 3 - 34-46 mmHg	PX-33		
Haddenham Custom GOLDPUNKT		Haddenham Custom GOLDPUNKT			
Class 2 - 23-32 mmHg	GP-22	Class 2 - 23-32 mmHg	GP-32	OPTIONS	
Class 3 - 34-46 mmHg	GP-23	Class 3 - 34-46 mmHg	GP-33		
Class 4 - over 49 mmHg	GP-24	Class 4 - over 49 mmHg	GP-34	PERTEX and GOLDPUNKT	CODE
THIGH HIGH + WAIST ATTACHMENT	7	PANTYHOSE		SILICONE GRIP TOP 2 Malleolous Pads Zip Velcro Fastening	CM-90 CM-91 CM-92 CM-93
Haddenham Custom PERTEX	CODE	Haddenham Custom PERTEX	CODE	VELCRO STRAP	CM-94
<b>Class 1</b> - 18-21 mmHg	PX-41	Class 1 - 18-21 mmHg	PX-51		
Class 2 - 23-32 mmHg	PX-42	Class 2 - 23-32 mmHg	PX-52		
Class 3 - 34-46 mmHg	PX-43	Class 3 - 34-46 mmHg	PX-53		
Haddenham Custom GOLDPUN Class 2 - 23-32 mmHg Class 3 - 34-46 mmHg Class 4 - over 49 mmHg	KT GP-42 GP-43 GP-44	Haddenham Custom GOLDPUN Class 2 - 23-32 mmHg Class 3 - 34-46 mmHg Class 4 - over 49 mmHg	KT GP-52 GP-53 GP-54		

Responsibility for the selection and application of products rests with the professional staff treating the patient.

Haddenham Healthcare does not take any responsibility for contraindicated use of these garments or for inaccurate measurements.